



Recombinant Human Palmitoyltransferase ZDHHC3 (ZDHHC3)

Product Code	CSB-CF882141HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NYG2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMLIPTHHFRNIERKPEYLQPEKCVPPYPGPVGTMWFI RDGCGIACAIVTWFL VLYAEF VVLVFM LIPSRDYVYSIINGIVFNLLAFLALASHCRAMLTDPGAVPKGNATKEFIE SLQL KPGQVVYKCPKCCSIKPDRAHHC SVCKRCIRKMDHHCPWVNNCVGENNQKY FVLFTMYIA LISLHALIMVGFHFLHCFEEDWTTYGLNREEMAETGISLHEKMQLNFSSTECS SFSPPT TVILLILLCFEGLLFLIFTSVMFGTQVHSICTDETGIEQLKKEERRWAKKTKWMN MKAVF GHPFSLGWASPFATPDQ GKADPYQYVV
Source	in vitro E.coli expression system
Target Names	ZDHHC3
Protein Names	Recommended name: Palmitoyltransferase ZDHHC3 EC= 2.3.1.- Alternative name(s): Protein DHHC1 Zinc finger DHHC domain-containing protein 3 Short name= DHHC-3 Zinc finger protein 373
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.