



Recombinant Human RING finger protein 175 (RNF175)

Product Code	CSB-CF019858HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N4F7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAGTAARKAAPVLEAPPQQEQLSHTKLSAEDTWNLQQERMYKMHRGHDS MHVEMILIFL CVLVIAQIVLVQWRQRHGRSYNLVTLLQMWVVPLYFTIKLYWWRFLSMWGMF SVITSYIL FRATRKPLSGRTPRLVYKWFLLIYKLSYAFGVVGYLAIMFTMCGFNLFFKIKAR DSMDFG IVSLFYGLYYGVMGRDFAEICSDYMASTIGFYSVSRLPTRSLSDNICAVCGQKII VELDE EGLIENTYQLSCNHVFHEFCIRGWCIVGKKQTCPYCKEKVDLKRMI SNP WERT HFLYGQI LDWLRYLVAWQPVVIGIVQGIISLGLE
Source	in vitro E.coli expression system
Target Names	RNF175
Protein Names	Recommended name: RING finger protein 175
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.