



Recombinant Human Germ cell-specific gene 1-like protein (GSG1L)

Product Code	CSB-CF740924HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6UXU4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKTSRRGRALLAVALNLLALLFATTAFLTTHWCQGTQRVPKPGCGQGGRANC PNSGANAT ANGTAAAPAAAAAATASGNGPPGGALYSWETGDDRFLFRNFHTGIWYSCEEE LSGLGEKC RSFIDLAPASEKGVWLWLSVVSEVLYILLLVVGFSLMCLELFHSSNVIDGLKLNAF AAVFT VLSGLLGMVAHMMYTQVFQVTVSLGPEDWRPHSWDYGWSFCLAWGSFTCC MAASVTTLNS YTKTVIEFRHKRKRKVFEEQGYREEPTFIDPEAIKYFRERMEKRDGSEEDFHLDGR HERYPAR HQPHMADSWPRSSAQEAPELNRQCWVLGHWW
Source	in vitro E.coli expression system
Target Names	GSG1L
Protein Names	Recommended name: Germ cell-specific gene 1-like protein Short name= GSG1-like protein
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.