



Recombinant *Saccharomyces cerevisiae* Palmitoyltransferase PFA3 (PFA3)

Product Code	CSB-CF337234SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42836
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MNDRLSLTSLFPRCLTTCLYIWTAYITLRIHQIPRWFLALTIVPTLAVALYTYK VIAR GPGSPLDFPDLLVHDLKAAENGLELPPEYMSKRCLTLKHDGRFRVCQVCHVW KPDRCHHC SSCDVCILKMDHHCPWFAECTGFRNQKFFIQFLMYTTLYAFLVLIYTCYELGTW FNLSGSF NRELIDFHLLGVALLAVAVFISVLAFTCFISIYQVCKNQTTIEVHGMRRYRRDLEIL NDSY GTNEHLENIFDLGSSMANWQDIMGTSWLEWILPIETFKYKKSHTKDEKGLYF NVRPQVQ DRLLSSRCLEDQLLRRVTPRPSLEADRASVEIIDAN
Source	in vitro E.coli expression system
Target Names	PFA3
Protein Names	Recommended name: Palmitoyltransferase PFA3 EC= 2.3.1.- Alternative name(s): Protein fatty acyltransferase 3
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.