



Recombinant Human Palmitoyltransferase ZDHHC15 (ZDHHC15)

Product Code	CSB-CF842720HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96MV8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRRGWKMALSGGLRCCRRVLSWVPVLVIVLVVLWSYYAYVFELCLVTVLSPA EKVIYLIL YHAIFVFFTWTYWKSIFTLPQQPNQKFHLSYTDKERYENEERPEVQKQMLVD MAKKLPVY TRTGSGAVRFCDRCHLIKPDRCCHHCSVCAMCVLKMDHHCPWVNNCIGFSNY KFFLQFLAY SVLYCLYIATTVFSYFIKYWRGELPSVRSKFHVLFLFVACMFFVSLVILFGYHC WLVSRL NKTTLAFCPTPVFTSGPEKNGFNLFGIKNIQQVFGDKKKFWLIPIGSSPGDGHS FPMRSM NESQNPLLANEETWEDNEDDNQDYPEGSSSLAVETET
Source	in vitro E.coli expression system
Target Names	ZDHHC15
Protein Names	Recommended name: Palmitoyltransferase ZDHHC15 EC= 2.3.1.- Alternative name(s): Zinc finger DHHC domain-containing protein 15 Short name= DHHC-15
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.