



Recombinant MscS family inner membrane protein Ynal (ynal)

Product Code	CSB-CF360085EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AEB6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MIAELFTNNALNLVIIFGSCAALILMSFWFRRGNRKRKRGFLFHAVQFLIYTIISAV GSI INYVIENYKLFITPGVIDFICTSLIAVILTIKLFLINQFEKQQIKKGRDITSARIMSR IIKITIIVVLVLLYGEHFGMSLSGLLTFGGIGGLAVGMAGKDILSNFFSGIMLYFDR PFS IGDWIRSPDRNIEGTVAEIGWRITKITTFDNRPLYVPNSLFSSISVENPGRMTNR RITTT IGLRYEDAAKVGVIVEAVREMLKNHPAIDQRQTLLVYFNQFADSSLNIMVYCFT KTTVWA EWLAAQQDVYLKIIDIVQSHGADFAFPSQTLYMDNITPPEQGR
Source	in vitro E.coli expression system
Target Names	ynal
Protein Names	Recommended name: MscS family inner membrane protein Ynal
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.