



Recombinant Human Dual oxidase maturation factor 1 (DUOXA1)

Product Code	CSB-CF007230HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1HG43
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MATLGHTFPFYAGPKPTFPMDDTLASIIMIFLTALATFIVILPGIRGKTRLFWLLRV VTS LFIGAAILAVNFSSEWSVGQVSTNTSYKAFSSEWISADIGLQVGLGGVNITLTGT PVQQL NETINYNEEFTWRLGENYAEYAKALEKGLPDPVLYLAEKFTPRSPCGLYRQY RLAGHYT SAMLWVAFLCWLLANVMLSMPLVYGGYMLLATGIFQLLALLFFSMATSLTSP CPLHLGA SVLHTHHGPAFWITLTTGLLCVLLGLAMAVAHRMQPHRLKAFFNQSVDEDPML EWSPEEG GLLSPRYRSMADSPKSQDIPLSEASSTKAYCKEAHPKDPDCAL
Source	in vitro E.coli expression system
Target Names	DUOXA1
Protein Names	Recommended name: Dual oxidase maturation factor 1 Alternative name(s): Dual oxidase activator 1 Numb-interacting protein
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.