



Recombinant *Oryza sativa* subsp. *indica* Photosystem Q (B) protein

Product Code	CSB-CF313923OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C433
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	TAILERRESTSLWGRFCNWITSTENRLYIGWFGVLMIPTLLTATSVFIIAFIAAPP VDID GIREPVSGSLLYGNNIISGAIPTSAAGLHFYPIWEAASVDEWLYNGGPYELIVL HLL GVACYMGREWELSFRLGMRPWIAVAYSAPVAAATAVFLIYPIGQGSFSDGMP LGISGTFN FMIVFQAEHNILMHPFHMLGVAGVFGGSLFSAMHGSLVTSSLIRETTENESAN EGYRFGQ EEETYNIVAAHGYFGRLIFQYASFNNSRSLHFFLAAWPVVGIWFTALGISTMAF NLNGFN FNQSVVDSQGRVINTWADIINRANLGMEVMHERNAHNFPDLA
Source	in vitro E.coli expression system
Target Names	psbA
Protein Names	Recommended name: Photosystem Q(B) protein EC= 1.10.3.9 Alternative name(s): 32 kDa thylakoid membrane protein Photosystem II protein D1
Expression Region	2-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.