



# Recombinant *Calycanthus floridus* var. *glaucus* Photosystem Q (B) protein (psbA)

<b>Product Code</b>	CSB-CF770311CBM
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q7YJY8
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Calycanthus floridus</i> var. <i>glaucus</i> (Eastern sweetshrub) ( <i>Calycanthus fertilis</i> var. <i>ferax</i> )
<b>Sequence</b>	MTAILERRESTSLWGRFCNWITSTENRLYIGWFGVLMIPTLLTATSVFIIAFIAAP PVDI DGIREPVSGSLLYGNIIISGAIPTSAAIGLHFYPIWEAASVDEWLYNGGPYELIV LHFL LGVACYMGREWELSFRLGMRPWIAVAYSAPVAAATAVFLIYPIGQGSFSDGM PLGISGTF NFMIVFQAEHNILMHPFHMLGVAGVFGGSLFSAMHGSLVTSSLIRETTENESA NAGYRFG QEEETYNIVAAHGYFGRLIFQYASFNNSRSLHFFLAAWPVVGIWFTALGISTMA FNLNGF NFNQSVVDSQGRVINTWADIINRANLGMEVMHERNAHNFPDLA
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	psbA
<b>Protein Names</b>	Recommended name: Photosystem Q(B) protein EC= 1.10.3.9 Alternative name(s): 32 kDa thylakoid membrane protein Photosystem II protein D1
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.