



Recombinant Danio rerio Protein YIPF3 (yipf3)

Product Code	CSB-CF772100DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q803Z2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSASQGSKNTNAEPWGGFDDNIIQGTGSAVIDMENMDDTSGSSFEDVGEMH QRMREEEEV TAEAAATEEDNGEYGEFLGMKGLKGQLGRQVADEVWQAGKRQASKAFNLYA NIDILRPYF DVEPIQVRNRLVESLIPVRMINFPQKVAGELYGPMMLVFTLVAILLHGMKTS GT VIREGT LMGTAIGTGFGYWLGVSFIYFLAYLCNAQITMLQMLSLLGYGLFGHCVVLFIT YNVHFH SLFYLLWMVIGGLSTLRMVAVLISRTVGQTPRLILCGSLAALHMLFLLYLHFAYH KMVEG ILDTLEGPNIPPIQRVARDVPVVASAVVNATVKSI AIVQSSQQL
Source	in vitro E.coli expression system
Target Names	yipf3
Protein Names	Recommended name: Protein YIPF3 Alternative name(s): YIP1 family member 3
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.