



Recombinant Human Protein FAM70A (FAM70A)

Product Code	CSB-CF708482HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5JRV8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MHQSLTQQRSSDMSLPDSMGAFNRRKRNSIYVTVLLIVSVLILTVGLAATTRT QNVTVG GYYPGVILGFGSFLGIIGSNLIENKRQMLVASIVFISFGVIAAFCCAIVDGVFAAR HIDL KPLYANRCHYVPKTSQKEAEEVISSSTKNSPSTRVMRNLTQAAREVNCPLHSR EFCTPRI RGNTCFCCDLYNCGNRVEITGGYYEYIDVSSCQDIIHLYHLLWSATILNIVGLFL GIITA AVLGGFKDMNPTLPALNCSVENTHPTVSYYAHPQVASYNYYHSPPHLPPYS AYDFQHSG VFPSSPPSGLSDEPQSASPSPSYMWSSSAPPRYSPPYPPFEKPPPYSP
Source	in vitro E.coli expression system
Target Names	TMEM255A
Protein Names	Recommended name: Protein FAM70A
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.