



Recombinant *Gossypium hirsutum* Photosystem II D2 protein (psbD)

Product Code	CSB-CF640742GHB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2L902
Product Type	Transmembrane Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Sequence	<p>MTIALGKFTKDENDLFDIMDDWLRDRFVFGWSGLLLFPCAYFALGGWFTG TTFVTSWY THGLASSYLEGCNFLTAAVSTPANSLSLAHSLLLLWGPEAQQDFTRWCQLGGLW TFVALHGA FGLIGFMLRQFELARSVQLRPYNIAIFSGPIAVFVSVFLIYPLGQSGWFFAPSF GVAAIF RFILFFQGFHNWTLNPFHMMGVAGVLGAALLCAIHGATVENTLFEDGDGANTF RAFNPTQ AEETYSMVTANRFWSQIFGVAFSNKRWLHFFMLFVPVTGLWMSALGVVGLAL NLRAYDFV SQEIRAAEDPEFETFYTKNILLNEGIRAWMAAQDQPHENLIFPEEVLPNGAL</p>
Source	in vitro E.coli expression system
Target Names	psbD
Protein Names	Recommended name: Photosystem II D2 protein Short name= PSII D2 protein EC= 1.10.3.9 Alternative name(s): Photosystem Q(A) protein
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.