



# Recombinant Secale cereale Photosystem II D2 protein (psbD)

<b>Product Code</b>	CSB-CF320794STH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P10803
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Secale cereale (Rye)
<b>Sequence</b>	<p>MTIALGRIPKEENDLFDTMDDWLRDRFVFGWSGLLLFPCAYFALGGWFTG  TTFVTSWY  THGLASSYLEGCNFLTAAVSTPANSLSLAHSLLLWGPEAQGDFTRWCQLGGLW  TFVALHGA  FALIGFMLRQFELARSVQLRPYNAISFSGPIAVFVSVFLIYPLGQSGWFFAPSFG  VAAIF  RFILFFQGFHNWTLNPFHMMGVAGVLGAALLCAIHGATVENTLFEDGDGANTF  RAFNPTQ  AEETYSMVTANRFWSQIFGVAFSNKRWLHFFMLFVPVTGLWMSAIGVVGLAL  NLRAYDFV  SQEIRAAEDPEFETFYTKNILLNEGIRAWMAAQDQPHENLIFPEEVLPRGNAL</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	psbD
<b>Protein Names</b>	Recommended name: Photosystem II D2 protein Short name= PSII D2 protein EC= 1.10.3.9 Alternative name(s): Photosystem Q(A) protein
<b>Expression Region</b>	1-353
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.