



Recombinant Human Rhomboid domain-containing protein 2 (RHBDD2)

Product Code	CSB-CF743985HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6NTF9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAASGPGCRSWCLCPEVPSATFFTALLSLLVSGPRLFLLQQPLAPSGLTLKSE ALRNWQV YRLVTYIFVYENPISLLCGAIIIWRFAGNFERTVGTVRHCFFTVIFAIFSAIIFLSFE AV SSLSKLGEVEDARGFTPVAFAMLGVTTVRSRMRRALVFGMVVPSVLVPWLLL GASWLIPQ TSFLSNVCGLSIGLAYGLTYCYSIDLSEVALKLDQTFPFSLMRRISVFKYVSGS SAERR AAQSRKLNVPVPGSYPTQSCHPHLSPSHPVSTQHASGQKLASWPSCTPGHM PTLPPYQPA SGLCYVQNHFQNPNTSSSVYPASAGTSLGIQPPTPVNSPGTVYSGALGTPGA AGSKESSR VPMP
Source	in vitro E.coli expression system
Target Names	RHBDD2
Protein Names	Recommended name: Rhomboid domain-containing protein 2
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.