



Recombinant Human Probable palmitoyltransferase ZDHHC20 (ZDHHC20)

Product Code	CSB-CF722800HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5W0Z9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPWTLWRCCQRVVGWVPVLFITFVVVWSYYAYVVELCVFTIFGNEENGKTV VYLVAFHL FFVMFVWSYWMTIFTSPASPSKEFYLSNSEKERYEKEFSQERQQEILRRAARA LPIYTTS ASKTIRYCEKCQLIKPDRAHHCSACDSCILKMDHHCPWVNNCVGFSNYKFFLL FLLYSLL YCLFVAATVLEYFIKFWTNELTDTRAKFHVLFVVSAMFFISVLSLFSYHCWLW GKNRT TIESFRAPTFSYGPDGNGFSLGCSKNWRQVFGDEKKYWLLPIFSSLGDGCSF PTRLVGMD PEQASVTNQNEYARSSGSNQPFPIKPLSESKNRLLDSESQWLENGAEEGIVKS GTNNHVT VAIEN
Source	in vitro E.coli expression system
Target Names	ZDHHC20
Protein Names	Recommended name: Probable palmitoyltransferase ZDHHC20 EC= 2.3.1.- Alternative name(s): Zinc finger DHHC domain-containing protein 20 Short name= DHHC-20
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.