



# Recombinant *Ashbya gossypii* Palmitoyltransferase ERF2 (ERF2)

<b>Product Code</b>	CSB-CF758696DOT
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q750R7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Sequence</b>	MRLHSRQASNPHRQYSAAQSLHSSSDSEHKEIPSTMGWAKRLMRWMVTVD QPHTFETSL KNYQSLAHVTNYIFFCGGRLRTVAKTKYLSVLVLMIAPIVLFVSVFETGYLWKH VAGAK PCVLCYYFWTLCFASFISTGATDPGTLPRNIHLAQLQDDYKLPLEYYSITLPS PVANA PVRLKYCTTCRIWRPPRASHCAVCDSCILSFDHHCDDLNNCIGQRNHRYFLAF LFSSVLS SIWLLTCCALKLRHAGSPSAAPVSLLLICYCAVSIWYPLLLAIYHLFLTGTQQTTH EYLK AVDSRNPIFHKVTHPERNPFVTGSCARNMLLLMCQPRGYDFLHTRSEHQAGD WRRFRLPI PHSFEKV
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	ERF2
<b>Protein Names</b>	Recommended name: Palmitoyltransferase ERF2 EC= 2.3.1.- Alternative name(s): DHHC cysteine-rich domain-containing protein ERF2 Ras protein acyltransferase
<b>Expression Region</b>	1-367
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.