



Recombinant Zea mays CASP-like protein 16

Product Code	CSB-CF471992ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6UBY6
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MASTPRTPAPERSPPPVPTPPPPLEDEPPPYLADGSPREEASFSSDGREGAP PKNPQLSP THHAAPRLVPPPSSPARQDGQEQEGSANKAAAATAEPAREPLRQMATGLAR PLSSQTSPA TTNSPTPSASPTPSSPAPVANNSKRSGQSTPKRAETKLPLSSPAVAVHFDPE EAVTSPL RLGKARLDQQQQQHAAGAAESGASVVPGVAAAVAERRELLALRLATA VLSLAAFS VIASARTSGWAGDYARHLQYRYAVAVNVIVFAYSVAQSLGKIRHLVSPRFTF RTMSSYY CSLFLDQVLAYLLMSASSAAASRNDLWVSRFGTDAFVRKITGALWLSFVAFLV LALNAVI SXANLFSMV
Source	in vitro E.coli expression system
Protein Names	Recommended name: CASP-like protein 16 Alternative name(s): Vegetative cell wall protein gp1
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.