



Recombinant Protein AdrA (adrA)

Product Code	CSB-CF359494EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AAP2
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MFPKIMNDENFFKAAAHGEEPPLTPQNEHQRSGLRFARRVRLPRAVGLAGM FLPIASTL VSHPPPGWWWLVLVGWAFVWPHLAWQIASRAVDPLSREIYNLKTDAVLAGM WVGVMGVNV LPSTAMLMIMCLNLMGAGGPRLFVAGLVLMVV SCLVTLELTGITVSFNSAPLE WWLSLPI IVIYPLLFGWVS YQTATKLAEHKRRLQVMSTRDGMTGVYNRRHWETMLRNEF DNCRRHNR DATLLIIDIDHFKSINDTWGHDVGDEAIVALTRQLQITLRGSDVIGRFGGDEFAVI MSGT PAESAITAMLRVHEGLNTRLRLPNTQVTLRISVGVAPLNPQM SHYREW LKSAD LALYKAK KAGRNRTEVAA
Source	in vitro E.coli expression system
Target Names	dgcC
Protein Names	Recommended name: Protein AdrA
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.