



# Recombinant Rabbit Gap junction alpha-1 protein (GJA1)

<b>Product Code</b>	CSB-CF754384RB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6TYA7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Sequence</b>	MGDWSALGKLLDKVQAYSTAGGKVVLSVLFIFRILLGTAVESAWGDEQSAFR CNTQQPGC ENVCYDKSFPISHVRFWVLQIIFVSVPTLLYLAHVFYVMRKEEKLNKKEEELKVA QTDGV NVEMHLKQIEIKKFKYGIEEHGKVKMRGLLRITYIISILFKSVFEVAFLLIQWYIY GFSL SAVYTCKRDPCPHQVDCFLSRPTEKTIFIFMLVSLVSLALNIIELFYVFFKGVK DRVK GKSDPYHATTGPLSPSKDCGSPKYAYFNGCSSPTAPLSPMSPPGYKLVTDGR NNSSCRNY NKQASEQNWANYSAEQNRMGQAGSTISNSHAQPFDFPDDNQNSKKLAAGHE LQPLAIVDQ RPSSRASSRASSRPRPDDLEI
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	GJA1
<b>Protein Names</b>	Recommended name: Gap junction alpha-1 protein Alternative name(s): Connexin-43 Short name= Cx43
<b>Expression Region</b>	1-382
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.