



Recombinant Escherichia coli Flagellar biosynthetic protein flhB (flhB)

Product Code	CSB-CF301195ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76299
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MSDESDDKTEAPTPHRLEKAREEGQIPRSRELTSLILLVGVSVIWFGGVSLAR RLSGML SAGLHFDHSIINDPNLILGQIILLIREAMLALLPLISGVVLVALISPVMLGGLVFSGK SL QPKFSKLNPLPGIKRMFSAQTGAELLKAILKTILVGSVTGFFLWHHWPQMMRL MAESPIT AMGNAMDLVGLCALLVVLGVIPMVGFDVFFQIFSHLKKLRMSRQDIRDEFKQS EGDPHVK GRIRMQRAAARRRMMADVPKADVIVNNPHTHYSVALQYDENKMSAPKVVAK GAGLVALRI REIGAENNVPTLEAPPLARALYRHAIEGQQIPGQLYAAVAEVLAWVWQLKRWR LAGGQRP VQPTHLPVPEALDFINEKPTHE
Source	in vitro E.coli expression system
Target Names	flhB
Protein Names	Recommended name: Flagellar biosynthetic protein flhB
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.