



Recombinant Drosophila melanogaster Innexin inx3 (inx3)

Product Code	CSB-CF896065DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VAS7
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MAVFGMVSAVSGFIKIRYLLDKAVIDNMVFRCHYRITTAIFTCCIIVTANNLIGDP ISC INDGAIPMHVINTFCWITYTYTIPGQQHRQIGTDVAGPGLGNEYGQEKRYHSYY QWVPFV LFFQGLMFYVPHWWKNMEDGKIRMITDGLRGMVSPDDYRRDRQDRILKYF VNSLNTHN GYSFAYFFCELLNFINVIVNIFMVDKFLGGAFMSYGTDVLKFSNMDQDKRFDP MIEIFPR LTKCTFHKFGPSGSVQKHDTLCVLALNILNEKIYIFLWFWFIILATISGVAVLYSLV VIM MPTTRETIKRSYRSAQRKEIAGLVRRLEIGDFLILHFLSQNLSTRSYSDMLQQL CGLLG ASRTPSAPSTLEMNRISHPIYPPVETFGGGKETET
Source	in vitro E.coli expression system
Target Names	Inx3
Protein Names	Recommended name: Innexin inx3 Short name= Innexin-3
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.