



Recombinant Uncharacterized protein ZK757.4 (ZK757.4)

Product Code	CSB-CF812830CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8I0G4
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSCRDCTRKGGCVYTTFWLVRFLPVVLVSLATGWGIYAYTYELCILSIDNWPQ RIIYLF FYALLILFYTSYLRTVYTKAWKPPQKYCIEGASKATYESVKDDERQLQLFLSDIA RERDL TLLVRGFDHGIRFCDKCCCIKPDRSHHCSMCEQCVLKFDHHCPWVNNCVNFG NYKYFILF LAYGFIFCIWIAATTLPSFIDFWRHEYDMNKKQYDSIDSVIQRNLKHLHTVLSNG RFPLV FLLFLSCMFSLSLFFYHLYLTAKNRTTVESFRAPMIDGKYAKDAFNHGIRAN YREIF GSHPLYWFLPVPSSIGDGCKFVMNDMTAMSAAAGNQVFVEMGNVPNGQSHA HPPIAQHHI NLYSEQSSSASEKEIKSVDDEITERSQVSTATTASPSHDSVTA
Source	in vitro E.coli expression system
Target Names	dhhc-4
Protein Names	Recommended name: Uncharacterized protein ZK757.4
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.