



Recombinant *Saccharomyces cerevisiae* Cytochrome b mRNA maturase bl2 (bl2)

Product Code	CSB-CF361172SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03873
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MAFRKSNVYLSLVNSYIIDSPQPSSINYWWNMGSLLGLCLVIQIVTGIFMAMHY SSNIEL AFSSVEHIMRDVHNGYILRYLHANGASFFFMVMFMHMAKGLYYGSYRSPRVT LWNVGVII FILTIATAFLGYCCVYGQMSHWGNMNIASNMFMNMKTIYMMMLMLLIYIFYTIM MRQMMK TKEYTMLIKSMDYINKNKYMINLNMTNKKDMNNNIGPLNMNLSIIYGSMLGDG HAEKRK GGKGRIVFQQEYCNINYLYLHSLLANLGYCNTNLPLIKTRLGKKGKIRQYLKF NTWTY DSFNMIYSEWYIKNMSGKGNKVIKSLDNYLTPLALAIWIMDDGCKLGKGLKFT TNCFS YKDVQYLTYLLHNKYNIKSTITKGNKENTQFVIYVWKESMPILTKIVSPYIIPSMK YKLG NYL
Source	in vitro E.coli expression system
Target Names	Bl2
Protein Names	Recommended name: Cytochrome b mRNA maturase bl2
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.