



Recombinant Mouse Voltage-dependent calcium channel gamma-8 subunit (Cacng8)

Product Code	CSB-CF823774MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VHW2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MESLKRWNEERGLWCEKGVQVLLTTIGAFSAFGLMTIAISTDYWLYTRALICNT TNLTAG DDGPPHRGGSGSSEKKDPGGLTHSGLWRICCLEGLKRGVCVKINHFPEDTDY DHDSA EYL LRVVRASSIFPILSAILLLLGGVCVAASRVYKSKRNIILGAGILFVAAGLSNIIGVIV YI SANAGEPGPKRDEEKKNHYSYGWSFYFGGLSFILAEVIGVLAVNIYIERSREAH CQSRSD LLKAGGGAGGSGGSGPSAILRLPSYRFRYRRRSRSSSRGSSEASPSRDASPG GPGGPGFA STDISMVTLSDPSKGSVAAGLASAGGGGSGAGVGAYGGAAGAAGGGGAGS ERDRGSSAG FLTLHNAFPKEAASGVTVTVTGPPAAPAPAPAPPAPAAPAGTLSKEAAAASNT NTLNRKT TPV
Source	in vitro E.coli expression system
Target Names	Cacng8
Protein Names	Recommended name: Voltage-dependent calcium channel gamma-8 subunit Alternative name(s): Neuronal voltage-gated calcium channel gamma-8 subunit Transmembrane AMPAR regulatory protein gamma-8 Short name= TARP gamma-8
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.