



Recombinant Mouse Pannexin-1 (Panx1)

Product Code	CSB-CF867490MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JIP4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAIAHLATEYVFSDFLAKEPTEPKFKGLRLELAVDKMVTICIAVGLPLLLISLAFQAQ EISI GTQISCFSPPSSFSWRQAAFVDSYCWAAVQQKSSLQSESGNLPLWLHKFFPYI LLLFAILL YLPALFWRFSAAPHLCSDLKFIMEELDKVYNRAIKAAKSARDLDRDGPDPGPPG VTENVGQ SLWEISESHFKYPIVEQYLKTKKNSSHLIMKYISCRLVTFVILLACIYLSYYFSL SLS DEFLCSIKSGVLKNDSTIPDRFQCKLIAVGIFQLLSLINLIVYALLIPVVVYTFPIF RQ KTDILKVYEILPTFDVLHFKSEGYNDLSLYNLFLEENISELKSYPKLVLENIKSN GQGI DPMLLLTNLGMKMDIIDGKIPTSLQTKGEDQGSQRVEFKDLDSSEAAANNGE KNSRQR LLNPSC
Source	in vitro E.coli expression system
Target Names	Panx1
Protein Names	Recommended name: Pannexin-1
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.