



Recombinant Escherichia coli UPF0053 inner membrane protein yfjD (yfjD)

Product Code	CSB-CF327748ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P37908
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MEHISTTTLLIILIMVVISAYFSGSETGMMTLNRYRLRHMAKQGNRS AKRVEKLL RKPD RLISLVLIGNNLVNILASALGTIVGMRLYGDAGVAIATGVLTFVVLVFAEVL PKTIA ALY PEKVAYPSSFLLAPLQILMMPLVWLLNAITRMLMRMMGIKTDIVVSGSLSKEEL RTIVHE SRSQISR RNQDMLLSVLDLEKMTVDDIMVPRSEIIGIDINDDWKSILRQLSHSPH GRIVL YRDSLDDAISMLRVREAWRLMSEKKEFTKETMLRAADEIYFVPEGTPLSTQLV KFQRNKK KVGLVVNEYGDIQGLVTVEDILEEIVGDFTTSMSP TLAEEVTPQNDG SVIIDGTA NVREI NKAFNWHLPEDDARTVNGVILEALEEIPVAGTRVRIGEYDIDILDVQDNMIKQVK VFPVK PLRESVAE
Source	in vitro E.coli expression system
Target Names	yfjD
Protein Names	Recommended name: UPF0053 inner membrane protein yfjD
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.