



# Recombinant Arabidopsis thaliana Reticulon-like protein B17 (RTNLB17)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-CF724961DOA  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | Q6DR04   |
| <b>Product Type</b>      | Transmembrane Protein  |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)   |
| <b>Sequence</b>          | MESTPPYHRSNTKSASRLQDSSNPPNLSLDLVLSSPNTPNPSSPVPLRDILLPLP<br>PSPLRK<br>SRTRLSDRLEMTSEDAMAVVRKRGKGGQKSLASPRNPRRSRRRSEAVE<br>EKEANLVIE<br>EIVKLPPRKRKTNGRPKKDKQSSAPPLCSSDLPNTCQSDLNLIGEIIISDLVMW<br>RDVAKS<br>TLWFGFGCLSFLSSCFAKGVNFSVFSAVSNLGLVLLCGSFLSNTLCQRKNEDT<br>KRAFHVS<br>EDDLRSARRVLPATNFFISKTSSELFSGEPSMTLKVTPFLLLGAEYGHILTLWRL<br>SAFGF<br>FLSFTIPKLYSCYTHQISQKVERVKTRIGEAWGVCSHKILAGSAVTAFWNLTSI<br>RTRIF<br>AVFIILVIFRYRRQNLQLTPEEVEPVENEQEEETLPQEEETVPQEEETVPQEEE<br>QTQPSE ERALVVVAET |
| <b>Source</b>            | in vitro E.coli expression system  |
| <b>Target Names</b>      | RTNLB17  |
| <b>Protein Names</b>     | Recommended name: Reticulon-like protein B17 Short name= AtRTNLB17   |
| <b>Expression Region</b> | 1-431  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |