



Recombinant Innexin-14 (inx-14)

Product Code	CSB-CF525326CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O62136
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MIWEIPIIGDFITKPFKAAKLYEFYDRLHLFTVYLLGFFVLLTGAKQHFGNPIDCM LPKQ HDDLKSWRDYIHNFLFYGTFRYDVSNGTSEFGSYTEDASVNYQWVPPFFA FQVCCFLL PFWCWAYMQKLIYIDMAFIVDYSKINSEKTFEKTKEKVDRIVNYMHDHFKFRR AHKMGY LSWITFNSAFPVLYSLTKLFFITNVIIQVNLVCKFLDVDSWTWGFDLLGKFIHPT PRAP EFSSFSDKQRFAAILTDGSYNRFQYFPILVGCEYQLQESVSNFVNHKAQCIIPM NVINEK IFIGLYFWLLVLTALSVIGTVKWILRIKSKKLNEMMIYKLIKKSLEREPFDSNIHDH RYN FVHKYLCADGILLIYFMMDTNGFLKTEEVIGALFKKYCSDAGLEPLQTAPTLTSS PRDLA YSETNYGFAAPQKI
Source	in vitro E.coli expression system
Target Names	inx-14
Protein Names	Recommended name: Innexin-14 Alternative name(s): Protein opu-14
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.