




Recombinant Ross river virus Structural polyprotein, partial

Product Code	CSB-CF357408RFM
Abbreviation	Recombinant Ross river virus Structural polyprotein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08491
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Ross river virus (strain T48) (RRV)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	YEHTATIPNVVGFPYKAHIERNGFSPMTLQLEVVETSWEPTLNLEYITCEYKTV VPSPFIKCCGTSECSSKEQPDYQCKVYTGVPFMWGGAYCFCDSENTQLSEA YVDRSDVCKHDHASAYKAHTASLKATIRISYG TINQTTEAFVNGEHAVNVGGS KFIFGP ISTAWSPFDNKIVVYKDDVYNQDFPPY GSGQPGRFGDIQSRTVESKD LYANTALKLSRPSPGVVHVPYTPTPSGFKYWLKEKGSSLNTKAPFGCKIKTNP VRAMDCAVGSIPVSM DIPDSAFTRVVDAPAVTDLSCQVVVCTHSSDFGGVATL SYKTDKPGKCAVHSHSNVATLQEATVDVKEDGKVTVHFSTASAPAFKVSVC DAKTTCTAACEPPKDHIVPYGASHNNQVFPDMSGTAMTWVQRLASGLGLGLAL IAVVVLVLVTCITMRR
Research Area	Immunology
Source	in vitro E.coli expression system
Target Names	N/A
Expression Region	817-1254aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	50.3 kDa
Protein Length	partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.