



Recombinant Danio rerio Glycerol-3-phosphate acyltransferase 3-like (agpat9l)

| | |
|--------------------------|--|
| Product Code | CSB-CF386475DIL |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A3KGT9 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Sequence | MEGYWAVLFPVLKVPWFSCVIVLIMLPAMFGISLGITETYMKLLIKTLEWATHRIQ RASRA EEILKESASNGLIQRDNSSLEQEIEELRRNRPKSADRGDFTLSDVLYFSRKGFE SIVEDD VTQRFTSEELVSWNLLTRTNNNFQYISLRLTVLWVVGWVRYCILLPLRITLTTI GLTWL VIGTTTTVGFLPNCRVKNWLSELVHLMCYRICARGLSATIHFNKQNRPKKGGIC VANHTS PIDVVILANDGCYAMVGQVHGGLMGVLQRAMERSCPHIWFERSEMRDRHLVT QRLKDHAVN AKTKLPILIFPEGTCINNTSVMFMFKKGSFEIGGTIYPVAIKYDPQFGDAFWNSSK YSIMS YLLRMMTSWAIVCNVWYLPPMTHEEGEDAVQFANRVKSTIAQQGGLVDLAW DGGLKRAKV KDSFKEQQQKKYSHMVVGEDSSD |
| Source | in vitro E.coli expression system |
| Target Names | agpat9l |
| Protein Names | Recommended name: Glycerol-3-phosphate acyltransferase 3-like EC= 2.3.1.15 Alternative name(s): 1-acylglycerol-3-phosphate O-acyltransferase 9-like Short name= 1-AGP acyltransferase 9-like Short name= 1-AGPAT 9-like Lysoph |
| Expression Region | 1-443 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |