



Recombinant *Candida albicans* Palmitoyltransferase PFA4 (PFA4)

Product Code	CSB-CF701624CZD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5AGV7
Product Type	Transmembrane Protein
Immunogen Species	<i>Candida albicans</i> (strain SC5314 / ATCC MYA-2876) (Yeast)
Sequence	<p>MAVQLKWPILGVIIIPCIIIFSLSYGSHYFILRHHLTMKQQLIYEFYVTMIWISYLLAI YT NPGRVPKKNYKPSLASSTRIEQTEDDSDGLGLESREDETLIREEPISGDRCEWIR YCKKCN NYKPPRSHHCKICQQCVLQMDHHCPWTLNLCVGNLPHFMRFLGWIIWGTG YLMIQLIK IINYENSNMPHYLFNKTELVAIIAITPLNFFVFASILVLFIRCLINICKGMTQIEIWE W ERLELQWSSKRLWRLIRFNYGRLHKGKPFPELNTWTNTTNNVNYNDNDDDG DEDVELINL ATNNNEDSTIVPQNFTIDDLIFPYNLGIWKNLVNALGYPYMWLIPFGKPKSNGY QPQISQ DYKQDDQLNLPWPPDGIRQKEIEINVLQQQGYQRDREEEEDDEELRSIRNYQEL RRRLDPR LNVQRSDFFINDMGEGLTDFGVDESD</p>
Source	in vitro <i>E.coli</i> expression system
Target Names	PFA4
Protein Names	Recommended name: Palmitoyltransferase PFA4 EC= 2.3.1.- Alternative name(s): Protein fatty acyltransferase 4
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.