



Recombinant UPF0053 inner membrane protein ytfL (ytfL)

Product Code	CSB-CF360055EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AE46
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	<p>MLNSILVILCLIAVSAFFSMSEISLAASRKIKLKLLADEGNINAQRVLNMQENPGM FFTV VQIGLNAVAILGGIVGDAAFSPAFLSLFSRYMSAELSEQLSFILSFSLVTGMFILF ADLT PKRIGMIAPEAVALRIINPMRFCLYVCTPLVWFFNGLANIIFRIFKLPVVRKDDIT SDDI YAVVEAGALAGVLRKQEHLEIENVFELESRTVPSSMTPRENVIVFDLHEDEQS LKNKVAE HPHSKFLVCNEDIDHIIGYVDSKDLLNRVLANQSLALNSGVQIRNTLIVPDTLTL EALE SFKTAGEDFAVIMNEYALVVGII TLNDVMTTLMGDLVGQGLEEQIVARDENSWL IDGGTP IDDMRVLDIDEFPQSGNYETIGGFMMFMLRKIPKRTDSVKFAGYKFEVVDIDN YRIDQL LVTRIDSKATALSPKLPDAKDKEESVA</p>
Source	in vitro E.coli expression system
Target Names	ytfL
Protein Names	Recommended name: UPF0053 inner membrane protein ytfL
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.