



Recombinant *Saccharomyces cerevisiae* Protein ECM7 (ECM7)

Product Code	CSB-CF312338SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06200
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MVMSRIRDITARPFFQNLTALEKVVQWLRLGTTLLIISFGLALTVGPLSSPRTLTM SRLDT YSADITTGLFTVLRESMEQSTSTEENNGVGLTTSELYILTAYTESQIKNVPQYIT VSLYG RCDSTYTMVEVFDSEGNMHSVKNSTTKSTCSSIGTDYLFDYREVLESGLDII DYAYNK IGSQQAESSAYTTYMRSLKHKKANVLHLLYAVISFQVCMLFFMIWYYYIKGRFM NALKER ALVHINSLLSLVVFIGGLISSISLAWVNYTIQSRINTELEAFGFSYHLGVTWFALL WCFA GLISVSCLAWSGLEWCISDNGTSYGGGIDDKFLGYQAGVFTDADLDDETSYSQ RYPQRQS TSGEALMRNSDTMATIRKTSVDLNSENDANTSLDHGNPTANISNGGKHEPF ATREEFE LQDIRFRSSNDSEESMQRVIKPSSALQF
Source	in vitro E.coli expression system
Target Names	ECM7
Protein Names	Recommended name: Protein ECM7 Alternative name(s): Extracellular mutant protein 7 Zinc-regulated gene 15 protein
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.