



Recombinant Mouse Protein tweety homolog 1 (Ttyh1)

Product Code	CSB-CF871973MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9D3A9
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MGAPPGYRPSAWVHLLHQLPRADFQLRPVPSGFAPRDQEYQQALLLVAALAG LGLGLSLI FIAVYLIRFCCCRPPEPHGAKSPPPGGCVTWSCIAALLVGCAGIGIGFYGNSE TSDGVS QLSSALLHANHTLSTIDDVLETVERLGEAVKTELTTLEEVLSVRMELVAATRG ARRQAE AAAQYLQGLAFWQGVSLSPVQVAEDVTFVEEYRWLAYVLLLLLVLLVCLFTLL GLAKQSK WLVVMTAMSLLVLSWGSMLAATAVGLSDFCSNPDTYVLNLTQEETGL SSDILSY FLCNQAVSNPFQQRLTSLQRALASIHSQLQGLEREAIPQFSAAQKPLLSLEETL NVTES FHQLVALLHCRSLHKDYGSALRGLCEDALEGLLFLMLFSLLSAGALATTLCSLP RAWALF PPSDDYDDTDDDDPFNPQESKRQVQWQSSI
Source	in vitro E.coli expression system
Target Names	Ttyh1
Protein Names	Recommended name: Protein tweety homolog 1 Short name= mTTY1
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.