



Recombinant Uncharacterized protein C07A9.8 (C07A9.8)

Product Code	CSB-CF327311CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34319
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MTNNNDPKCVVIEEPDDVPVKPPKIFDFTDWPFEIPDIFKRKELSYNNYDLAT SKSLMI VRMIFKWRGSLWQAVYKELIVWICAYSLVSVIYRFALTRSQKDIFERFGEYCDA RMGYLP LNFVLGFFCNIIRRWLKLTYSLGNIDNIALFVSAYVRGTDDRARQIRRNIIRYCVI SQC LVFRDIHVGVRRRFPTLEAVAQAGIMLPHELEKFNSIKSRYQKYWVSFNWALE LLNVAKT EKSIDGDNARNIAAQEISKFRSALTTVSMYDWVPIPLMYPQLVNMAVHTYFFLC IFTRQF FISADAHNKTEVDLYIPFMTIIEFIFYMGWLKVAMELLNPFGEADDDFDCNLLID RNLAI GLTSVDDAYDQLPEVKPDVFTGGSVKPLDSDDTRSLKYHFGSAAQMEEISYLK KEENKMI AAGKKPNKLLWVKSRRRKFETSATQPSFPIP
Source	in vitro E.coli expression system
Target Names	best-5
Protein Names	Recommended name: Uncharacterized protein C07A9.8
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.