



Recombinant Ectopic P granules protein 3 (epg-3)

Product Code	CSB-CF893184CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9XWU8
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MAKKQKKSTEKSERTVEFKEPPKPANSEERLKPAGRGMKPSPSQNTLNRME RETIVFWRR PHIVIPYALMEIAHLAVELFFKILAHKTVLLLTAISIGLAVYGYHAPGAHQEHVQTI EKH ILWWSWWVLLGVLSSIGLGSGLHTFLIYLGPHIAAVTMAAYECQSLDFPQPPYP ESIQCP STKSSIAVTFWQIVAKVRVESLLWGAGTALGELPPYFMARAARISGQEPDDEE YREFLEL MNADKESDADQKLSIVERAKSWVEHNIHRLGFPGILLFASIPNPLFDLAGITCGH FLVPF WSFFGATLIGKALVKMHVQMGFVILAFSDHHAENFVKILEKIPAVGPYIRQPISD LLEKQ RKALHKTPGEHSEQSTSYLAWGLSLMVTFMILFFFLSIVNSLAKDYHKRLWER KRRQNKD LIDEENQSFEEEEEEAVTPPSSCPLLLSDGFEGVVVKK
Source	in vitro E.coli expression system
Target Names	epg-3
Protein Names	Recommended name: Ectopic P granules protein 3
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.