



Recombinant Uncharacterized protein C03F11.1 (C03F11.1)

Product Code	CSB-CF606124CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q11122
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MKHEQRKHSDFWKRRKISQSAMVTQCESNGSVTSHNTSSAFQRNNSRYGVP IDSTAVKQV YRVRSERLTNRVRITDRSLYLALFGVILMLVESEITAEKFYGVSKTHWISQSLRV GVTCS TIALLYHIILYHLNDIVLELVDCGADDWRVVVTTTERVIQFCIEFICCGICPLPGSGE MKW TFIEPSLHKDGPSEERTVQTRNVDVPVDVILSCFMLCRSYLFARFMVLHASKQF QDASTRT LAALNRIQVNFVSVIKTSLDQQPVFLTTFTFIFWIIMSWMFVQCERYGFSGKNP QSILY SNSLWFIAITFMLNGYGDIVPQTNAGRFAIFVGVVGAVISSILIAVISRNILLSQG QRN VNNFMYDSKLAREHKEAAARVLQHTWHIHKCLQGSDDGNRRLRTYQRKFLKA IHKFRSVK SEMREFSENNSQNNPQMTRLVTDMTSMQRLNVDQEMRMQIEVVKFLK
Source	in vitro E.coli expression system
Target Names	kcnl-3
Protein Names	Recommended name: Uncharacterized protein C03F11.1
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.