



Recombinant Rat Acetylcholine receptor subunit epsilon (Chrne)

Product Code	CSB-CF005400RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09660
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>KNEELSLYHHLFDNYDPECRPVRRPEDTVTITLKVTLTNLISLNEKEETLTTSVW IGIEW QDYRLNFSKDDFAGVEILRVPSEHVWLPEIVLENNIDGQFGVAYDCNVLVYEG GSVSWLP PAIYRSTCAVEVTYFPFDWQNC SLIFRSQTYNAEEVELIFAVDDDGNAINKIDID TAAFT ENGEWAIDYCPGMIRHYEGGSTEDPGETDVIYTLIIRRKPLFYVINIIVPCVLISGL VLL AYFLPAQAGGQKCTVSINVLLAQTVFLFLIAQKIPETSLSVPLLGRYLIFVMVAT LIVM NCVIVLNVSLRTPPTTHATSPRLRQILLELLPRLGLSPPPEDPGAASPARRASSV GILLR AEELILKKPRSELVFEGQRHRHGTWTAALCQNLGAAPEVRCCVDAVNFVA ESTRDQEA TGEELSDWVRMGKALDNVCFWAALVLFVSVGSTLIFLGGYFNQVPDLPYPPCIQ P</p>
Source	in vitro E.coli expression system
Target Names	Chrne
Protein Names	Recommended name: Acetylcholine receptor subunit epsilon
Expression Region	21-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>