



Recombinant Arabidopsis thaliana Probable S-acyltransferase At3g56930 (At3g56930)

Product Code	CSB-CF878668DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M1K5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	<p>MAWNETKLKRLYQVWRGSNKFLCGGRLIFGPDASSLYLSTILILGPAVMFFVK MYTKMAD PRTKNPNLCIPILCVSWILTILDIFFLLMTSSRDPGIVPRSFRRPETDDAPDSTTP SMEW VSGRTPNIRIPRVKDVTVNGHTVKVKFCDTCLLYRPPRASHCSICNNCVQRFD HHCPWVG QCIGVRNYRFFFMFISTSTTLCIYVFAFSWLNIFQRHMDEKISIWKAISKDVLSDI LIVY CFITVWFVGGTLIFHSYLICTNQTTYENFRYRYDKKENPYNKGILGNIWEIFLSKI PPSM NKFRSFVKEEDYMMMMVETPTS NLGESLVSSKEKIDIEMGGGRIVDESGKSYS LPEILRN LNYEDLEDDCEEDDLKAKDHHHHHHHQHQHNEGIIPFPDFFFTNEIGSNKDER NGEESGG SSSDGENTGKRVRVSEDEEEKVEGYERNWSTDKGMNINAGSEDGASSPVST SPMLRK</p>
Source	in vitro E.coli expression system
Target Names	PAT04
Protein Names	<p>Recommended name: Probable S-acyltransferase At3g56930 EC= 2.3.1.- Alternative name(s): Probable palmitoyltransferase At3g56930 Zinc finger DHHC domain-containing protein At3g56930</p>
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>