



Recombinant Human Neuronal acetylcholine receptor subunit beta-2 (CHRNA2)

Product Code	CSB-CF005396HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17787
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	TDTEERLVEHLLDPSRYNKLIRPATNGSELVTVQLMVSLAQLISVHEREQIMTT NVWLTQ EWEDYRLTWKPEEFDNMKKVRLPSKHIWLPDVVLYNNADGMYEVSFYSNAVV SYDGSIFW LPPAIYKSACKIEVKHFPFDQQNCTMKFRSWTYDRTEIDLVLKSEVASLDDFTP SGEWDI VALPGRRNENPDDSTYVDITYDFIIRRKPLFYTINLIIPCVLITSLAILVFYLPSCDG EK MTLCSVLLALTVFLLLISKIVPPTSLDVPLVGKYLMTMLVTFIVTSVCVLNVH HRS PTTHMAPWVKVVFLEKLPALLFMQQPRHHCARQRLRLRRRQREREGAGALF FREAPGAD SCTCFVNRASVQGLAGAFGAEPAPVAGPGRSGEPCGCGLREAVDGVRFIAD HMRSEDDQ SVSEDWKYVAMVIDRLFLWIFVFCVFGTIGMFLQPLFQNYTTTTFLHSDHSAP SSK
Source	in vitro E.coli expression system
Target Names	CHRNA2
Protein Names	Recommended name: Neuronal acetylcholine receptor subunit beta-2
Expression Region	26-502
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.