



# Recombinant *Schizosaccharomyces pombe* Probable acyl-CoA desaturase (SPCC1281.06c)

<b>Product Code</b>	CSB-CF528462SXV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94523
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Sequence</b>	MTAPSATAFSSATTQPTTEGNASMRKRTIPVVPSVPERKWDPKAPKHIQEQP WTMQNWWR HLNWLHCMLIFGLPMIAIYGVFTTPLQTKTLIFAIYYAYSGLGITAGYHRLWSHR AYKA KKPLEYFLAAGGAAAFEGSIRWWSRDHRAHHRYSRDTDKDPYINVKKGFWYAH VGWMILQN PRRIGRSDVSDLNSDPFVMFNHRHFLPIASFMAFIFPSLFCGLLWGDYRGGYF YAGVCRL VVFVHHATFCVNSLAHLIGSQPFDDTNSARNHFITALVTLGEGNHNYHHAFPND YRNGLRW YEYDPTKIFIYIASLFLAYNLNTPDNEIQKGIVQQKQKVLDRWRACLNWGIPL EQLPV MEFEDFLEQSKTRPLVLINGVVHDMTGFEHPGGQGLLRSFAGKDATAAFNGG VYDHTNGA HNLLSTYRVAVVRGGMEVEVWKSGAGAQLPMKDTQGQKIVRVGEQITRIQPPI EAAAAN
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SPCC1281.06c
<b>Protein Names</b>	Recommended name: Probable acyl-CoA desaturase EC= 1.14.19.1 Alternative name(s): Delta(9)-desaturase Short name= Delta-9 desaturase Fatty acid desaturase Stearoyl-CoA desaturase
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.