



Recombinant Sphingomyelin synthase-related 1 (F53B1.2)

Product Code	CSB-CF633621CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q20696
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MLDNRPIPADPNEWRCEDVGNWLKKIGMAKYADLIAMKHKVDGKCLLALTD DLKDPPVS INCLGDIKKILFAIEFLSQKVVVEIGNSGVHHRSTPNGNGPSLKNSKDGLLVEYNE QNHLS ISGEDVYTTTRRAEIVEDEETLLDTLAKSSDGTSTVQLISREEIIRQVERPDTYFK SVAK LLIAFAYSSLSFLMTSFVMVLVHDRVDPDTKTYPLPDIVLDNVPHIPWAFDMCE TIGLVL AVVWFTVLFHNRVIVARRMFSLGTVFLLRCFTMLITSLVPGIHLQCEARP NTTMQE KLHKAFHIWSNLGMSLHGVRSCGDYMFSGHTTVITMISHFITEYTPADWTGLH TFTWVLN CFAIFLILAAHEHYSIDVFIAFYISSRMFLYYHAYAYNHAGITATDYRMRTWFPLG WFFE YGSQGVENEFSLPINIRIPRRVFFAKSEEPKITPKSDSSRKRSSVVAAKQNGN SKNHTK KHN
Source	in vitro E.coli expression system
Target Names	F53B1.2
Protein Names	Recommended name: Sphingomyelin synthase-related 1
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.