



Recombinant Bacillus subtilis Uncharacterized ABC transporter permease yclI (yclI)

Product Code	CSB-CF309413BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P94412
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	<p>MNFIKRAFWNMKAKKGKTLQLFVFTVICVFVLSGLAIQSAAQKSSELARQELG GSVTLQ VDRQKQMEKQQDSGEKRSFESTPIKVSDANKLAALDHSVKSINYTTTSASANAG NFDAIESS SSSDSSSSSSSSNAKNSQGGGQGGPQMVQADLSIEGVISTALVDDFSDGDSK ITDGRAIT KSDVGKKVTVINETLAEENDLSVGDSITIESATDEDTTVCLKIVGIYKTTSSGDD QAQNF SFLNPYNKLYTPYTATAALKGDDYKNTIDSAVYYMDDAKNMDTFVKAACKTSID FDTYTL NTNDQLYQQMVGPIENVASFKNVVYLVSVAGAVILGLIVMMSIRERKYEMGV LMAIGEK RWKLIQQLTEILIVAVIAIGLASVTGNLVANQLGNQLLSQQISSSTDSTQTASG QMPGG GGGMGGKMFHSSSNVDVIDSLNVAVSMNDMLILGGIGILIAIIATLLPSISVLRLL HPKT ILTKQE</p>
Source	in vitro E.coli expression system
Target Names	yclI
Protein Names	Recommended name: Uncharacterized ABC transporter permease yclI
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>