



Recombinant Chicken Gap junction alpha-3 protein (GJA3)

Product Code	CSB-CF009446CH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29415
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	GDWSFLGRLLENAQEHSTVIGK VWLTVLFIFRILVLGAAAEVWGDEQSDFTC NTQQPGCENVCYDKAFPIHSHIRFWVLQIIFVSTPTLIYLGHVLHIVRMEEKRKEK EEELKKRGSVKDNNYPGAATSGGGSGGGNFKDPPIKMGKEKLPIRDERRI RMGGALLRTYIFNIIFKTLFEVGFIVGQYFLYGFELKPVYQCSRPPCPHTVDCFI SRPTEKTIFIIFMLVVASVSLLLNMLEIYHLGWKCLKQGMTSQYSLEMPVTTLTP VMVTGESKPVSLPPPAPPVVVTTTAPAPVLPDTRAVTPLLAPVTMAPYAAAA PRTRPPSNTASMASYPVAPPV PENRHRAVTPTPVSTPVTIPTPIPTPTPAIINYF NSKSNALAAEQNWVNMAAEQQGKAPSSSAGSSTPSSVRHPLPEQEEPLEQL LPLPAGPPITTTNSGSSTSLSGASGSKWDVEGEEELAEERPISATCTTVMEMHE PPLLVDTRRLSRASKSSSSRARSDDLAV
Source	in vitro E.coli expression system
Target Names	GJA3
Protein Names	Recommended name: Gap junction alpha-3 protein Alternative name(s): Connexin-56 Short name= Cx56
Expression Region	2-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.