



Recombinant Arabidopsis thaliana ALBINO3-like protein 3, mitochondrial (ALB3L3)

Product Code	CSB-CF605888DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0WUC5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	<p>STPSDLDELTRFRDDSIAGVGSNDHGLEFDDSIASLGSNHGLEFGDSIGSL GSNGYGL EFGDLSQDLSNYDYLTQPVISLLDSYHDITGLPWWVVIATSTVAFRTALLPILILQ RKQT KRISQFLPKLPHFWPPQGSGRSVLDQLKLFKRDIGCPSFLWVPAYFSIQIS CFFLWI TSIRRMSLDHHPGFDSSGGALWFQNLTEIPNGLYGPLFPFLIAGLHYTNTQITFT ASSVHK VDKFTELAKAYKTFNLNLLTCALYFLSFQMPQGSLLYWATNLSFSIAQQSILNHP VVSAKL GLQANDSVQKEAGNPILTININEGKLTDPSSKGRLLSVHNLTPKELVALSAKYLS GGHKDK SIPLRLALEKDPEYLQAMIILGQALYQKDQFAEAAKCLEQAASKLLDTSPTVEVE EVDLL IVASQWAGVSNIRQGKTSEGITHLERNVANMKEPDDPKSKAHYLDALVLYSSAIF NEGRRE EAAKYLRVAVYDPSFSELLKQCEEDDTIPTSSSSNSTSKTS</p>
Source	in vitro E.coli expression system
Target Names	ALB3L3
Protein Names	Recommended name: ALBINO3-like protein 3, mitochondrial
Expression Region	45-566
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.