



Recombinant Innexin unc-7 (unc-7)

Product Code	CSB-CF312700CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q03412
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MLGSSSNPEPPLLSRIIGVPPPPPPRAPTTALVLRVPTVDSPSKKKQQPDTRNK YQETAL RDKKTRTPLEKARHLDNLPSYQAQKLLDGSHQLRIDSHHVGSAGHGAGQGHHG HKKEFGPA MILYYLASAFRALYPRLDDDFVDKLNYYYYTTTILASFALLVSAKQYVGFPIQCWV PATFT DAMEQYTENYCWVQNTYWVPMQEDIPREIYSRRNRQIGYYQWVPFILAIEALL FYVPCIL WRGLLYWHSGINLQGLVQMACDARLMDSEIKTRTVYTMARHMQDEVQLTNID RQGHRSRSC FSNLQLGANCGRHCGCYVTMLYIGIKVLYSANVLLQFFLLNHLLGSNDLAYGFS LLKDLM HAIEWEQTGMFPRVTLCDFEVRVLGNIHRHTVQCVLMINMFNEKIFLFLWFWF LTCGIVT VCNTMYWILIMFIPSQGMSFVRKYLRVLPDHPAKPIADDVTLRKFTNNFLRKDG VFMLRM ISTHAGELMSSELILALWQDFNNVDRSPTQFWD AEHGQGTID
Source	in vitro E.coli expression system
Target Names	unc-7
Protein Names	Recommended name: Innexin unc-7 Alternative name(s): Uncoordinated protein 7
Expression Region	1-522
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.