



Recombinant *Drosophila virilis* Lethal (2)neighbour of Tid protein (I (2)not)

Product Code	CSB-CF630509DMP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q24332
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Sequence	<p>MAPPRTNRRPVTRRKKPSNLLLLLKEKYVNWNFYKYLAFEPAAALPIVSVLIVLAE AVINV LVIQRVPYTEIDWKAYMQECEGFLNGTTNYSLLRGDTGPLVYPAAFVYIYSGLY YLTGQG TNVRLAQYIFACIYLLQMCLVRLRYTKSRKVPPYVLVLSAFTSYRIHSIYVLRFLN DPVA ILLLYAALNLFDRWTLGSLCYSLAVGVKMNILLFAPALLLFYLANLGVLRTLVQ LTIC AVLQLFIGAPFLRTHPMEYLRGSFDLGRIFEHKWTVNYRFLSKELFEQREFHLV LLALHL LLLLCFAKPTWTFRRSYLRLREVQQQIVPQILHKNREQQNQAKIQAKKSKQLVE QQQEE DEQLSAEQKSFLKAFERSLQQSTGQKSLEPHLNQKQEYDIYFEPCTQLALLP FFLCNFI GVACARSLHYQFYIWFHSLPYLVWSTPYSGLGVRYLILGIIYCWNTYPSTMES SIALHL CHLILLVGVARHVYRILRLNATVKQQQRQEQLQNGQPKTSKKSQ</p>
Source	in vitro E.coli expression system
Target Names	I(2)not
Protein Names	Recommended name: Lethal(2)neighbour of Tid protein Alternative name(s): NOT58
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.