



Recombinant Mouse Vang-like protein 1 (Vangl1)

Product Code	CSB-CF773945MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q80Z96
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MDTESTYSGYSYYSSHSKSHRQGERTREERHKSPRNKDGRGSEKSVTIQAPA GEPLLAND SARTGAEEVQDDNWGETTTAITGTSEHSISQEDIARISKDMEDSVGLDCKRYL GLTVASF LGLLVFLTPIAFILLPQILWREELKPCGAICEGLLISVSFKLLILLIGTWALFFRKQR AD VPRVVFVFRALLLVLIFLVVSYWLFYGVRIILDSRDQNYKDQVAVSLVDALLFIH YLAI VLLELRQLQPMFTLQVVRSTDGESRFYSLGHLSIQRAALVVLENYYKDFTIYNP NLLTAS KFRAAKHMAGLKVYNVDGSPNATGQSRAMIAAAARRRDSSHNELYYEEAEH ERRVKKRR ARLVVAVEEAFIHIQRLQAEQQKSPGEVMDPREAAQAIFPSMARALQKYLRT TRQQHYH SMESILQHAFICITNSMTPKAFLERYLSAGPTLQYDKDRWLSTQWRLISEEAVT NGLRDG IVFVLKCLDFSLVNVKIPFIVLSEEFIDPKSHKFLRLQSETSV
Source	in vitro E.coli expression system
Target Names	Vangl1
Protein Names	Recommended name: Vang-like protein 1 Alternative name(s): Loop-tail protein 2 Van Gogh-like protein 1
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.