



# Recombinant *Oryza sativa* subsp. japonica Probable mannan synthase 9 (CSLA9)

<b>Product Code</b>	CSB-CF714993OFG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q67VS7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Sequence</b>	<p>           MAAAGAVLPEQIAAMWEQVKAPVVVPLLRLSVAACLAMSVMLFVEKVMYSVV            LVGVHLFG            RRPDRRYRCDPIVAAGADNDDPELADANAAPFMVLIQIPMYNEREVYKLSIGA            ACGLSWP            SDRVIVQVLDDSTDPVIKEMVQVECKRWESKGVRIKYEIRDNRVGYKAGALRE            GMKHGYV            RDCDYVAIFDADFQDPDFLARTIPFLVHNPDIQALVQARWKFVNANECLMTRM            QEMSLDY            HFKVEQEVGSSTHAFFGFNGTAGVWRISAMNEAGGWKDRTTVEDMDLAVRA            GLKGWKFVY            LGDLMVKSELPSTFKAFRYQQHRWSCGPANLFRKMLVEIATNKKVTLWKKIYV            IYNFFLV            RKIIGHIVTFVYCLVVPATVLIPEVEIPRWGYVYLPISIVTILNSIGTPRSLHLLIFW            VL            FENVMSLHRTKATLIGLLETGRVNEWVTEKLGDAKLLKLPKAFRRPRMRIG            DRVNALE LGFSAYLSFCGCYDIAYGKGYSSLFLFLQSITFFIIGVGYVGTIVPH         </p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CSLA9
<b>Protein Names</b>	Recommended name: Probable mannan synthase 9 EC= 2.4.1.- Alternative name(s): Cellulose synthase-like protein A9 OsCslA9
<b>Expression Region</b>	1-527
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>