



Recombinant Danio rerio Protein tweety homolog 2-like (ttyh2l)

Product Code	CSB-CF757516DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6NUZ2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MASSRQDYIAPWWTYWLHNPHLNFNFQTVDNTEFKPEDASYQQSLVFLACVS AVALGLCL LLLSVYLTCLCCCRREEDEEVKRPDTCCVTWAAVITGLVICSAVGVGFGNSE TNDGVYQ LTYSIYNANHTLGGIGSMVGNLSLGMQVGLKEHLERLDEIFSPRGDYTQTLRF MQQMADN IIQQLTAMPETRKVQVDLAAIADRTAFIEYYRWLTYLLLLLIDLVICLLACLALAKQ SRW LLTVIMVCGMLTLIMSWASLGAGTATAVGTSDFCVSPDKYIVNQTKGTLSSDIV HYYLYC SQSLPNPFQQSLTIFQRSLTMMQIQVQGLLQFAVVVFPTAEKDLLGIQRLNNS EFNLHQ LTALLDCRSLHKDYLEALLGVCYDGVGELLYLCLFSLAACAFCALLCAVPRAW MLIAIR DRDYDDIDEEDPFNPTRRFNAYNPSRAQVHSFCSYSSSMGSQTSLQPPLQAT SSTPEYMN QSMLFERNPRYENVPLIGSGSPPPSYSPSMRTTYLSMTDAQIRHFGTDFQV</p>
Source	in vitro E.coli expression system
Target Names	ttyh2l
Protein Names	Recommended name: Protein tweety homolog 2-like
Expression Region	1-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.